

Prothrombin Time (PT) Reagent Lot Number Change

Audience: Internal and/or External Clients
Issuing Division: Division of Hematological Pathology
Enquiry: Jacqueline King (416) 480-4299 or Rita Selby (416) 480-4048
E-mail: Jacqueline.King@sunnybrook.ca or Rita.Selby@sunnybrook.ca

Effective April 8, 2019

The Division of Hematological Pathology will be changing the lot number of the Thromborel S (Prothrombin Time) reagent and as a result the current reference range will change.

The new reference ranges are listed below for the PT reference range:

Test Name	New Age Range	New Reference Range	Previous Age Range	Previous Reference Range	Units
PT	>3 months	10.3-12.9	>3 months	10.7-13.3	seconds
	0-3 months	10.3-17.9	0- 3 months	10.7-18.5	seconds

Please refer to patient reports for details. If you have any comments, concerns, or require further information, do not hesitate to contact Dr. Rita Selby.

The information in this Bulletin is current as of the date indicated above. However, it may be superseded by subsequent communications. Please check to make sure you are reviewing the most recent communique."