

## Prothrombin Time (PT) Reagent Lot Number Change

**Audience:** Internal and/or External Clients  
**Issuing Division:** Division of Hematological Pathology  
**Enquiry:** Jacqueline King (416) 480-4299 or Rita Selby (416) 480-4048  
**E-mail:** [Jacqueline.King@sunnybrook.ca](mailto:Jacqueline.King@sunnybrook.ca) or [Rita.Selby@sunnybrook.ca](mailto:Rita.Selby@sunnybrook.ca)

### Effective June 18, 2020

The Division of Hematological Pathology will be changing the lot number of the Thromborel S (Prothrombin Time) reagent and as a result the current reference range will change.

The new reference ranges are listed below for the PT reference range:

Test Name	New Age Range	New Reference Range	Previous Age Range	Previous Reference Range	Units
PT	>3 months	10.0-12.5	>3 months	10.3-12.9	seconds
	0-3 months	10.0-17.1	0- 3 months	10.3-17.9	seconds

Please refer to patient reports for details. If you have any comments, concerns, or require further information, do not hesitate to contact Dr. Rita Selby at 416-480-4048 or email: [Rita.Selby@sunnybrook.ca](mailto:Rita.Selby@sunnybrook.ca).